




 POLARIS

PROFESSIONAL SERVICES
CORPORATION
10000 POLARIS DRIVE
MINNETONKA, MN 55343
TEL: 952.892.2000
WWW.POLARIS.COM



POLARIS

Dental Implant System 

GC, Dental Implant System

POLARIS

Finds the Way...

Introduction	2
Standards	4
Philosophy behind Polaris	6
Macro-mechanics	6
Micro-mechanics	8
Critical Exercises	11
Portfolio	14
Closure	14
Appendices & Auxiliaries	18
Surgical Plan	27

Introducing the Implant Project

This company started a road plan of "Manufacturing the Group of Dental Implant Products" in 2016 and localized the knowledge of manufacturing this product with the help of the educated youth and scientific cities and transferred technology from the top companies in the world. The company applied the invaluable plan in the field of health industry in 2020 through employing the world's labor and the most advanced equipment and number one European raw materials and by supplying the specialized installations and infrastructure equipment, implementing international systems and technical standards, and through effective cooperation with the accredited international universities and international reference laboratories. The company's portfolio includes a broad spectrum of the group of equipment and materials in the field of dentistry and dental prosthetics. It will be improved by introducing this new product and the company will get one step closer to its main goal, which is supplying the dental technology and commodity requirements to its fullest extent.

Dental implant surgery has emerged as a highly effective and reliable solution for replacing the roots of missing natural teeth, with success rates exceeding 95% over 10 years in clinical studies. This proven track record has led to a significant increase in the adoption of implantology worldwide, making it the gold standard for tooth replacement.

The success of an implant system depends on several critical factors, including the design of the implant's external surface, the type of connection, the coating technology, the shape of prosthetic components, and the position of the vaginal lip. Here, we provide a detailed overview of the Polaris B&P portfolio, a comprehensive and innovative solution designed to meet the highest standards in dental implantology, ensuring optimal outcomes for both clinicians and patients.



Standards





NARROW

Connections:

Morse III[®]

Hex: 1.7

MLA Thread

Other Involving
Structural Features:

- Open 2° axial bevel with tapered chamfer below the shoulder
- Flat lower bevel
- chamfered chamfered profile

Material:

Material:

- 316L stainless steel
- Titanium (commercially pure titanium)

Hybrid Design:

- Special type

Cutting Edge:

- Cutting edge
- chamfering

Our portfolio includes implants in Narrow and Regular lines, with the Regular line also featuring short implants for smaller clinical applications. Additionally, a Wide line will soon be introduced, maintaining the same connection type as the Regular line for consistency and ease of use. The Narrow line features a 1.7 mm hex, while the Regular line is designed with a 2.0 mm hex, ensuring precision and adaptability for a wide range of clinical scenarios.

REGULAR

Connection:
Double Morse T²
Hex: 2.5
M2 Thread

Coronal

Coronal (e.g. **Coronal**)
 Coronal (e.g. **Coronal**)
 Coronal (e.g. **Coronal**)
 Coronal (e.g. **Coronal**)


Straight Body

Tapered Apex

Tapered Apex (e.g. **Tapered Apex**)
 Tapered Apex (e.g. **Tapered Apex**)



To complement our implant systems, we offer a comprehensive range of prosthetic components designed to meet diverse clinical needs. Our portfolio includes standard prosthetic parts, digital solutions for streamlined workflows, and customized elements tailored to individual patient anatomy for optimal aesthetics and functionality.



Enhanced tissue integration: The porous structure promotes tissue integration through increased cell adhesion and bone formation around the implant.

Surface roughness (R_a - 0.8 micrometers): The treatment creates both macro and micro roughness on the implant surface. This roughness increases the surface area and energy, which helps in better bone-implant contact (BIC). The porous surface allows for better mechanical interlocking with the bone, enhancing the primary stability of the implant.

Reduced healing time: The porous structure allows for a shorter healing period. The optimized surface topography allows bone cells to colonize the implant more quickly and accelerates the healing process and reduces surface coating.

Increased durability/hydrophobicity: The porous structure improves the wettability of the implant, which enhances the contact between the implant and the biological environment, further supporting osseointegration.

5.2.2 SURFACE

Surface roughness and porosity



What we are looking for from the manufacturer's viewpoint is the fine surface as different contact area will create proper morphology for seal and preventing smaller imperfections such as 500µm indicates the effect of oxidation and higher imperfections shows the pure effect of thermal treatment. Most of the voids found are smaller in size or they become the secondary imperfections in the surface oxide.

State-of-the-Art Manufacturing Infrastructure

The foundation of high quality dental implant manufacturing is a structured infrastructure and process design process. At the facility, we have established modern manufacturing lines equipped with cutting edge technology to ensure the highest standards of accuracy, consistency, and reliability throughout the production process. Within these systems, for greater efficiency, we highly automate measurement systems for quality control, and a certified machine for optimal surface texture. Additionally, we employ a semi-robotic surface treatment process to enhance implant bio-compatibility and performance. To ensure the utmost hygiene and sterility, our manufacturing takes place in a certified cleanroom environment, ensuring that every implant meets stringent international standards. The robust infrastructure enables a modular system that is not only flexible within production but allows to meet the diverse needs of clinicians and patients worldwide.



Critical Properties of Dental Implants



The success of dental implants hinges on their fundamental aspects: physical properties, chemical properties, and cellular response. The physical properties, such as mechanical strength, shape, roughness, hydrophobicity, and mechanical strength, ensure the implant can withstand functional loads and integrate seamlessly with surrounding bone.

From the chemical composition viewpoint, the implant material must exhibit high biocompatibility and inherent resistance to prevent cellular reactions in the biological environment.

Most importantly, the cellular response to the implant—be it toxicity, surface characteristics, and material composition—plays a pivotal role in osseointegration, ensuring stable and long-lasting integration with the patients. By optimizing these properties using the results of SEM, cell adhesion, and cell differentiation tests, we create implants that not only perform reliably but also promote natural healing and bone formation.





Packaging Labeling

- Double Sealing
- Including isolation Kapoids & Filter Layer
- Or you prefer Sterilization



-  Lot number
-  Expiry date
-  Certificate of origin included
-  No other supply date
-  For single use only
-  Good sterility control
-  Temperature controlled
-  Keep away from heat
-  Keep away from water
-  Performance



- Always open the package by cutting the square from the top.

BLP Implant



Wisely Chosen, Simply Done...

Series Length	Series Diameter								
	Series (mm 1.5)			Series (mm 2)				Series (mm 3)	
	Ø1.5	Ø1.6	Ø1.8	Ø2.0	Ø2.2	Ø2.5	Ø2.8	Ø3.0	Ø3.2
1.50						☒	☒	☒	☒
1.80		☒	☒	☒	☒	☒	☒	☒	☒
1.50-2	☒	☒	☒	☒	☒	☒	☒	☒	
1.80-2	☒	☒	☒	☒	☒	☒	☒	☒	
1.50-3	☒	☒	☒	☒	☒	☒	☒		







Prosthetics & Auxiliaries



Wisely Chosen, Simply Done...

Prosthesis Diagram



Cover Screen

Following steps:
 1. Cover screen
 2. Heat in the oven



Applicable for Superfine Resin
 1. 100°C/212°F
 2. 10-15 minutes

Heating Abutment

Following steps:
 1. Heat in the oven



Impression Coping Pickup



Impression Coping Transferor



Type		diameter	length	Ref.
Impression				
Pickup	Short	4.0	15.0	OPH-4015
	Long	4.0	19.0	OPH-4019
Transferor	Short	4.0	8.0	CTP-4008
	Long	4.0	16.0	CTP-4016
Transferor				
Pickup	Short	4.0	15.0	OPH-4015B
	Long	4.0	19.0	OPH-4019B
Transferor	Short	4.0	8.0	CTP-4008B
	Long	4.0	16.0	CTP-4016B

Straight Abutment Cementable



Abutments für Implantate von 3,30 mm

Artikelnummer	Ø	Ø2	Ø3	Ø4	Ø5
02010	3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30
02011	3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30
02012	3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30
02013	3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30
02014	3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30
02015	3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30

Abutments für Implantate von 4,10 mm

Artikelnummer	Ø	Ø2	Ø3	Ø4	Ø5	Ø6
02020	4,10	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30
02021	4,10	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30
02022	4,10	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30
02023	4,10	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30
02024	4,10	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30
02025	4,10	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30
02026	4,10	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30	Ø3,30/3,30



Angled Abutment Connectable



Abutment Height Options (mm)



Orange Series

	10.0	12.0	15.0	18.0	22.0
Ø-17° 08 02 15	ØK1000108	ØK1000112	ØK1000115	ØK1000118	ØK1000122
Ø-17° 08 02 15	ØK1000108	ØK1000112	ØK1000115	ØK1000118	ØK1000122



Green Series

	10.0	12.0	15.0	18.0	22.0
Ø-17° 08 02 15	ØK1000108	ØK1000112	ØK1000115	ØK1000118	ØK1000122
Ø-17° 08 02 15	ØK1000108	ØK1000112	ØK1000115	ØK1000118	ØK1000122
Ø-17° 08 02 15	ØK1000108	ØK1000112	ØK1000115	ØK1000118	ØK1000122
Ø-17° 08 02 15	ØK1000108	ØK1000112	ØK1000115	ØK1000118	ØK1000122
Ø-17° 08 02 15	ØK1000108	ØK1000112	ØK1000115	ØK1000118	ØK1000122
Ø-17° 08 02 15	ØK1000108	ØK1000112	ØK1000115	ØK1000118	ØK1000122



Orange Series
Abutment Connectable



Green Series
Abutment Connectable

UCLA Abutment

10000



10000



Temporary Abutment

10000



10000



Ball Abutment

10000




10000



Prosthetic Diagram

Ball Abductor Orientation



Fine Millod Abrasives

Model	System	Dimensions	Part Code
A-200	Reverse	50	2010001-001
A-201	Reverse	50	2010001-002
A-202	Reverse	50	2010001-003

Model	System	Dimensions	Part Code
A-203	Angular	50	2010001-004
A-204	Angular	50	2010001-005
A-205	Angular	50	2010001-006



Small Series

Large Body

Model	System	Part Code
A-206	Reverse	2010001-007
A-207	Angular	2010001-008



Small Series

1/2 Size

Model	System	Part Code
A-208	Reverse	2010001-009
A-209	Angular	2010001-010



Small Series



Surgical Plan

simple & friendly

PMMA® Surgical Kit



BLT Drill


Hardening with G3 Rockwell

Sharp drill with 3 Flute

Coated with TiN

Life length up to 200 surgeries





Typical Procedure

Lightweight concrete aggregate (LCA) or
lightweight aggregate (LWA)



—————●—————
Compressive Test

—————○—————
Tensile Test



Lightweight concrete aggregate (LCA) or
lightweight aggregate (LWA)



—————●—————
Compressive Test

—————○—————
Tensile Test




Abbildung 10: Anwendungsbeispiel 10 (10):
 Abstr. Abstr. (10) (10)




Abbildung 11: Anwendungsbeispiel 11 (11):
 Abstr. Abstr. (11) (11)

