



Auto MIX (Fully Automatic)



QUALITY & ADVANCED TECHNOLOGY
PROVIDE FULL FACILITIES
AS YOU DESIRED

Auto MIX (Fully Automatic)

This device is used in dental laboratories in order to mix materials (plasters used in dentistry) and also to uniform and remove air bubble.



Features

- Microprocessor pressure control system
- Internal oil-less vacuum pump
- Widespread timer & program adjustment Time
- Digital Indicator in order to display vacuum surface with adjustment capability
- The possibility to change time parameter & operation
- Vacuum termination valve to enhance function safety
- Shut down alarm system
- The possibility to modify and change functional programs
- Easy consumable parts replacement
- Easy usage and installation
- Light-weighted & portable



Technical Data

- Dimensions: 28 cm (width)*28 cm (depth)*48 cm (height)
- Weight: 17 Kg
- Maximum power input: 85 W
- Operating voltage: 210-240 V-50/60 Hz
- Micromotor Angular velocity: 400 RPM
- Maximum amount of vacuum: 870 mBar
- Container volume: 250 and 500 ml

Address: Humboldt Street 12 - 040723
 Hilden - Germany
 Tel: + 49 (0) 21033373500
 Fax: + 49 (0) 21033373550
 Email: info@polarisgmbh.de
www.polarisgmbh.de

24
 Months Warranty

10
 Years Supporting of
 Spare Parts