



Dual Dental Laboratory Work Bench



QUALITY & ADVANCED TECHNOLOGY
PROVIDE FULL FACILITIES
AS YOU DESIRED



Dual Dental Laboratory Work Bench



Features

- Built-in micromotor system, controllable foot pedal, adjustable speed, clockwise & counter-clockwise rotation
- Strong air nozzle to clean the dusts
- Suction system to collect suspended particles & dust
- Particle collection bags with filtering valves
- Automatic on & off suction system controlled by micromotor
- Standard fluorescent lights with adjustable light intensity
- Metallic stands with ability to adjust height
- Fast & easy selection and replacement of table stands according to working environment
- Heat-resistant, anti-scratch & anti-bacterial surface (made from composite stone or MDF)
- Shelf (console) for purpose of placing the necessary instruments
- Two elbow support to enhance ergonomics
- Protective shield in order to enhance safety
- Three separated drawers and two underneath drawer with plastic container (tray)
- Internal gas, air & water connections
- Benzene flame with internal gas connections
- Electrical Switchgear with four outgoing socket
- The possibility to install 40000 RPM micromotor with elbow control
- Programmable intelligent system and screen display to control and display parameters (micromotor rotational speed, light intensity) in a modular package

Technical Data

- Dimensions: 60 cm (width)*110 cm (depth)*86-100 cm (height)
- Power input of dust evacuation system: 1400 W
- Operating voltage: 210-240 V-50/60 Hz
- Micromotor Angular velocity: 30000-45000 RPM
- Air pressure: 1-6 Bar

Address: Humboldt Street 12 - D40723
Hilden - Germany
Tel: + 49 (0) 21033373500
Fax: + 49 (0) 21033373550
Email: info@polarisgmbh.de
www.polarisgmbh.de

24
Months Warranty

10
Years Supporting of
Spare Parts