



QUALITY & ADVANCED TECHNOLOGY
PROVIDE FULL FACILITIES
AS YOU DESIRED

Flexoplast:

Flexoplast is an advanced system for injecting thermoplastic material with ideal usage in dentistry to build Acrylic removable prosthetics (denture). Easy, precise and elegant prosthetics operation is result of Flexoplast using.



! Features

- Equipped with an advanced processor to automatically control program
- Equipped with Program termination alarm
- Alarm system for faulty user operation
- Precise thermal sensor to measure the muffle temperature
- Indicator to display muffle temperature
- Made of High quality cast iron to facilitate casting in cartridge injection & ejection
- High precision performance
- User friendly

● Cartridge Features

- High polishability and inexistence of saliva absorption
- With high biocompatibility to adapt to oral tissues
- Fabricated in four colors: Transparent, Intense Pink, Soft Pink, Dark Pink
- Fabricated in Large, Medium and Small sizes

⚙️ Technical Data

- Dimensions: 20 cm (width)×30 cm (depth)×55 cm (height)
- Weight: 30 Kg
- Power input: 800 W
- Operating voltage: 210-240 V-50/60 Hz

Address: Humboldt Street 12 - D40723
Hilden - Germany
Tel: + 49 (0) 21033373500
Fax: + 49 (0) 21033373550
Email: info@polarisgmbh.de
www.polarisgmbh.de

24
Months Warranty

10
Years Supporting of
Spare Parts