



Pin-hole drilling unit

Pindex Q1



QUALITY & ADVANCED TECHNOLOGY
PROVIDE FULL FACILITIES
AS YOU DESIRED

Pindex:

Pin-hole drilling unit with laser pointing system for precise model drilling required by the removable pinstechnique adjustable drilling depth.



Features

- Adjustable drilling depth
- Stable, Resistant, Balanced Design Three tips for High accuracy
- User Friendly
- Intelligent timer for optimized consumption
- Contains Powerful and noiseless engine
- Accurate and reliable point laser system
- Safety operation by using mechanical and electronic protective equipment



Technical Data

- Dimensions: 25 cm (width)*26 cm (depth)*38 cm (height)
- Technical data: Maximum power input: 150 W to 2800 RPM
- 3mm shank pliers for drill
- Net weight: 18 Kg
- Power supply: 220V-60/50 Hz



Address: Humboldt Street 12 - 040723
 Hilden - Germany
 Tel: + 49 (0) 21033373500
 Fax: + 49 (0) 21033373550
 Email: info@polarisgmbh.de
www.polarisgmbh.de

24
 Months Warranty

10
 Years Supporting of
 Spare Parts