



# Steam Cleaner



**QUALITY** & ADVANCED TECHNOLOGY  
PROVIDE FULL FACILITIES  
AS YOU DESIRED



# Steam Cleaner

## Steam Cleaner:

Steam Cleaning is a prevalent and effective method to clean dental models, prosthesis and tools, removing casting operation and residual wax in dental and dentistry. Steam Cleaner of Polaris with easy and safety operation has been manufactured to fulfil maximum requirements of dental and dentistry.



### Features

- Eliminate pollution of different surfaces for metallic and ceramic (Zirconium) frames
- Cleaning different surfaces of prostheses
- Steam adjusting output according to the sort of consumption
- Three tips for multiple consumption cases
- Water amount electronic controller
- Water amount alarm system
- Portable and fixed nozzle
- Fully metallic framework for strengthening and enhancement of product life purposes
- Steam storage device for standby status
- Safety valves to reach maximum operational security.
- Foot pedal



### Technical Data

- Dimensions: 30 cm (width)\*46 cm (depth)\*40 cm (height)
- Maximum output steam pressure: 6 Bar
- Net weight: 18 Kg
- Maximum power input: 2200 W
- Power supply: 210-240 V-50/60 Hz
- Water storage volume: 4 liters

Address: Humboldt Street 12 - 040723  
Hilden - Germany  
Tel: + 49 (0) 21033373500  
Fax: + 49 (0) 21033373550  
Email: [info@polarisgmbh.de](mailto:info@polarisgmbh.de)  
[www.polarisgmbh.de](http://www.polarisgmbh.de)

**24**  
Months Warranty

**10**  
Years Supporting of  
Spare Parts