



# Sunny Therm - 1



**QUALITY** & ADVANCED TECHNOLOGY  
PROVIDE FULL FACILITIES  
AS YOU DESIRED

## Sunny Therm - 1

The application of this product is to pre-warm and prepare casting cylinders for alloy frames of chrome-cobalt and porcelain. Also, this furnace is used in pre-warming partial frame casting cylinders.



### Features

Performing heating and dewaxing for chrome cobalt, nickel, chromium prosthesis cylinders and other precious and non-precious alloys

Containing 10 functional programs to change all parameters  
Programmable temperature with three distinct stages (time, temperature and temperature increase rate)

Ideal case with high capacity for setting 4 cylinders with weight of 180 g

Weekly schedule (100 hours)

LCD for displaying all parameters

Microcontroller for precise & real-time control

Intelligent alarm system to warn software/hardware errors

Precise thermal sensor (thermocouple type K)

Quick installation & ease of operation

Using the best ceramic fiber thermal insulators

Using the best thermal elements with high durability

### Technical Data

Dimensions: 38 cm (width) \* 38 cm (depth) \* 48 cm (height)

Chamber Dimensions: 16 cm (width) x 16 cm (depth) x 10 cm (height)

Chamber capacity: 4 cylinders

Maximum temperature: 1150° C (2100° F)

Maximum power consumption: 1800W

Weight: 25 kg

Power supply: 210-240 V AC -50/60 Hz

Address: Humboldt Street 12 - 040723  
Hilden - Germany  
Tel: + 49 (0) 21033373500  
Fax: + 49 (0) 21033373550  
Email: [info@polarisgmbh.de](mailto:info@polarisgmbh.de)  
[www.polarisgmbh.de](http://www.polarisgmbh.de)

**24**  
Months Warranty

**10**  
Years Supporting of  
Spare Parts