



# Vibrator



**QUALITY** & ADVANCED TECHNOLOGY  
PROVIDE FULL FACILITIES  
AS YOU DESIRED

## Vibrator

The application of Vibrator in industry is preparing uniform powder. In dental field, Vibrator is used to create uniform surfaces in plasters, all kinds of refractory and casting materials and air bubble removal. Vibration adjustment and easy operation are considered as sensible characteristics of this product.



### Features

- Uniform long life vibration system
- Vibration level adjustment system
- An ideal product for elimination of bubbles from plaster
- Applicable in other industries for a similar application for other powders
- Light-weighted & portable along with easy Usage



### Technical Data

- Dimensions: 16 cm (width)\*24 cm (depth)\*12 cm (height)
- Net weight: 7 Kg
- Maximum power input: 100 W
- Power supply: 210-240 V-50 Hz

Address: Humboldt Street 12 - D40723  
Hilden - Germany  
Tel: + 49 (0) 21033373500  
Fax: + 49 (0) 21033373550  
Email: [info@polarisgmbh.de](mailto:info@polarisgmbh.de)  
[www.polarisgmbh.de](http://www.polarisgmbh.de)

**24**  
Months Warranty

**10**  
Years Supporting of  
Spare Parts